

HOLDING micropipettes

Holding micropipets are used to hold the oocyte or blastocys

Microtech standard Holding micropipettes are prepared from borosilicate glass tubing (O.D. 1 mm, I.D. 0,75 mm), total length 5,5 cm, polished opening, length of arm 0,9 mm, bending angle 20°- 40°. They may be ordered in 3 different sizes small O.D. 80 µm, medium O.D. 100 µm and large O.D. 120 µm. Small and medium have polished opening an I.D. of 15-20 µm, large size has an I.D. 25 µm. The micropipettes are made strait or with bending angle.



Individually packed

For single use only

20 pcs in box

Identity colour

GREEN	Article number	code	O.D.	angle	Polished opening
small					
	001-80-20	001	80	20	15
	001-80-30	001	80	30	15
	001-80-35	001	80	35	15
	001-80-40	001	80	40	15
	001-80	001	80	strait	15
medium					
	001-100-20	001	100	20	20
	001-100-30	001	100	30	20
	001-100-35	001	100	35	20
	001-100-40	001	100	40	20
	001-100	001	100	strait	20
large					
	001-120-20	001	120	20	25
	001-120-30	001	120	30	25
	001-120-35	001	120	35	25
	001-120-40	001	120	40	25
	001-120	001	120	strait	25

Bending angle



Strait

